


# Medika 140

## Vaccine Refrigerator

**fiocchetti**  
SCIENTIFIC REFRIGERATORS

Made in  
Italy



### SPECIFICATION

Dimensions (WxDxH)	520x635x835 mm
Temperature range (alarmed)	+2°C to +8°
Capacity	140 Litre
Reversible glass door	Standard
Door lock with key	Standard
Internal lighting	Standard
Audible Door Alarm	Standard
Alarm Contacts	Standard
N° of shelves	2
Forced air ventilated cooling	Standard
Automatic evaporation of defrost condensate	Standard
Automatic defrosting	Standard
Power/Currenty	140 Watts/1.4 Amp
Refrigerant gas	R134A
Temperature probe port (15mm)	Standard
Aluminum pull out drawers	Optional

BELL TECHNOLOGY LTD  
10c Maurice Rd, Penrose, Auckland  
PO Box 26353, Epsom, Auckland

Tel: (64) 9 5251875  
Fax: (64) 9 525 1874  
Email: info@belltechnology.co.nz

**Bell**  
TECHNOLOGY LTD