

# Medika 250

Vaccine Refrigerator

**fiocchetti**  
SCIENTIFIC REFRIGERATORS

Made in  
Italy



## SPECIFICATION

Dimensions (WxDxH)	600x655x1635 mm
Temperature range	+2°C to +8°
Capacity	250 Litre
Reversible glass door	Standard
Door lock with key	Standard
Internal LED lighting	Standard
Audible Door Alarm	Standard
Alarm Contacts	Standard
N° of shelves	4
Forced air ventilated cooling	Standard
Automatic evaporation of defrost condensate	Standard
Automatic defrosting	Standard
Castors	Standard
Temperature probe port (15mm)	Standard
Aluminum pull out drawers	Optional

BELL TECHNOLOGY LTD  
10c Maurice Rd, Penrose, Auckland  
PO Box 26353, Epsom, Auckland

Tel: (64) 9 5251875  
Fax: (64) 9 525 1874  
Email: [info@belltechnology.co.nz](mailto:info@belltechnology.co.nz)

**Bell**  
TECHNOLOGY LTD