

EMOTECA 140



Characteristics

Structure	White or Stainless steel 18/10 int.+ext.
Insulation thickness	40 (polyur. mm)
Doors number	1
Door type	Glass
Refrigeration type	Ventilated
Refrigerant type	R134A
Temperature range	+4°
Temperature control/adjustment	ECT-F Control® : with Hi/Low Temp. alarm, door ajar, ON/OFF w/password
Defrosting type	Auto
Voltage	230V/1/50Hz
Castors	standard fitted
Capacity (LT)	140
Safety key locking	standard fitted
Internal lighting	standard fitted
DMLP Digital Monitor (Alarm System)	standard fitted
Drawers nr.	3
Storage capacity (blood bags nr., 450 ml)	63
Width (cm)	52
Depth (cm)	64
Height (cm)	100
Net weight (kg)	72
Equipment in accordance with	Directive 93/42/EEC annex V, article 3 and following modifications
Automatic water evaporation Probe for product simulation Ext. terminal board for RS485/dry contacts Available also with fully insulated door	

Accessories

Drawer perspex front (H=83 mm)	DMLP PT100 Extension module + nr.1 PT100 probe + 4-20mA output
PT100 probe	Frequency 60Hz
WEB LIGHT SERVER (Ethernet + GSM)	Voltage 115V
GSM Communicator CF2008 module	DMLP Printer
DMLP Remote panel	